



SPECIFICATION SHEET

HIP & JOINT Gel / Paste for Dogs & Cats

Glucosamine & Chondroitin & MSM

- Hip & Joint Inflammation
- Muscle Problems

DOG & CAT Label Ingredients: Water, Soy Bean Oil, Stearic Acid, Lecithin, Molasses, Poultry Flavor, Glucosamine, Sodium Chondroitin Sulfate, Methylsulfonylmethane, Yucca Schidigera Extract, Sea Kelp Extract, Xanthan Gum, Sodium Carboxymethyl Cellulose, Ascorbic Acid, Lactic Acid, Sodium Benzoate, Potassium Sorbate. (Ingredients may change without notice.)



Analytical Parameters (Actives)	per 6 grams....	per 100 grams
Glucosamine (Shellfish Sourced).....	350mg.....	5830 mg
Sodium Chondroitin Sulfate (Bovine Sourced).....	50mg.....	830mg
MSM.....	250mg.....	4165mg
Vitamin C Ascorbic Acid.....	20mg.....	330mg
Yucca Schidigera Extract.....	15mg.....	250mg
Sea Kelp Extract.....	1mg.....	16mg



Analytical Parameters (Inactives)	per 6 grams....	per 100 grams
Crude Protein.....0.5%.....	30mg.....	500 mg
Crude Fat.....22%.....	1320mg.....	22g
Crude Fiber.....0.5%.....	30mg.....	500mg
Moisture.....74%.....	4400mg.....	74g
Proprietary Blend: Vegetable Soy Bean Oil, Stearic Acid, Lecithin, Molasses, Poultry Flavor, Salmon Oil, Xanthan Gum, Sodium Carboxymethyl Cellulose, Ascorbic Acid, Lactic Acid, Sodium Benzoate, Potassium Sorbate, Water.....quantity sufficient to 100% pH 3.5 to 5		

(Specifications may change without notice)

3 OUNCE TUBE	5 OUNCE TUBE
White tube	White tube
White dispensing cap	White dispensing cap
Tube size 1.5 x 6 inches	Tube size 1.5 x 7 inches
Carton qty. 12	Carton qty. 12
carton weight 3.25 pounds	carton weight 4.75 pounds
carton size 6 x 6 x 6 inches	carton size 6 x 6 x 8 inches

We researched the optimal support for helping our four-legged friends to stay active. Here it is. Tasty and easy to give.

Directions for Use

Apply 1/4 teaspoon to your pet's paw or bowl once or twice daily. Amount can be adjusted according to size of animal.

Heavy Metals Parameters

Arsenic (As).....	NMT 12 ppm
Cadmium (Cd).....	NMT 10 ppm
Mercury (HG).....	NMT 10 ppm
Lead (Pb).....	NMT 0.27 ppm

Microbiological Parameters

Total Plate Count.....	NMT 1,000 CFU/g
Yeast and Molds.....	NMT 100 CFU/g
Coliform.....	Absent in 10g
Salmonella spp.....	Absent in 10g
E. coli.....	Absent in 10g
S. aureus.....	Absent in 10g

(NMT no more than; NLT no less than)